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
The significance of viscoelasticity in foodstuff rheology

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The rheological properties of foodstuffs are of considerable importance. Often, however, only simple, standard methods are used to characterise the substances, e.g. the flow and tension test. Detailed assessments of the substances are only, however, possible by carrying out vibration determinations. The only way the food rheologist can find answers to questions relating to storage stability, retrogradation, regeneration and sedimentation is through dynamic determinations. In this paper, the connection between non-destructive determinations and practice-related problems is demonstrated.

1 Fundamental principles of viscoelasticity

The most interesting phenomena are exhibited by viscoelastic fluids which behave like elastic solids and viscous liquids. The appearance of an elastic component means that some of the energy is used not for shear but for a reversible elastic deformation. The cause of this is the orientation or deformation and destruction of the aggregates or of a structure through shear. A viscoelastic fluid can be clearly characterised by means of the vibration test, in which the substance is subjected to a harmonically oscillating shear deformation. The deformation γ is given by equation (1), where γ_0 is the deformation amplitude and $\omega = 2\pi r$ is the angular frequency.

As a response of the system one obtains a vibrating shear stress τ at the same angular frequency ω .

This simple mathematical relationship only, however, applies in the linear-viscoelastic region [1]. The sample can be elastically deformed up to a maximum deformation $\gamma_{0,max}$, as in the case of Hooke's Law. Above this maximum deformation the viscosity will be dependent upon the shear rate, which can only be described by non-linear differential equations which cannot be solved or, if so, then only approximately.

By forming the differential coefficient (equation 2) one obtains equation (3) from (1).

For ideal viscous solutions the shear stress τ is directly proportional to the shear rate $\dot{\gamma}$, which is equivalent to a phase displacement of 90° (see equation 3).

$$\gamma = \gamma_0 \sin(\omega t) \quad (1)$$

$$\frac{d\gamma}{dt} = \dot{\gamma} \quad (2)$$

$$\dot{\gamma} = \omega\gamma_0 \cos(\omega t) \quad (3)$$

oder/or

$$\dot{\gamma} = \omega\gamma_0 \sin(\omega t + \delta), \delta = 90^\circ$$

A purely elastic fluid is linearly connected with the deformation, via Hooke's Law and produces the phase

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displacement $\delta = 0^\circ$ (see equation 1). A viscoelastic solution therefore has a phase displacement of between 0° and 90° .

If one defines the viscosity and the shear modulus as a complex quantity [1, 2, 3], one obtains, taking equations (1) and (2) into account, the equations (4) and (5) with

The storage modulus G' is a measure for the elastically stored energy, whilst the loss modulus G'' represents the energy which is dissipated through viscous flow. From these two quantities it is possible to calculate the loss factor $\tan \delta$ according to equation (6).

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2 Applications

2.1 Sedimentation

Every disperse system is exposed to shear forces. The „particles“ sink to the bottom or apparently rise to the surface, depending on the density of the dispersed substance. The force acting on a particle having a volume V and density ρ can be calculated according to equation (7), where g is the acceleration due to gravity and ρ_L the density of the dispersion medium.

In order to keep the particles in suspension therefore, a force which acts against the sedimentation force has

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was born in 1957 and completed his physics studies between 1978 and 1983, graduating between 1984 and 1987. In 1986 he founded Suck Wissenschaftliche Geräte Entwicklung, the purpose being to correlate macroscopic flow behaviour with microscopic aggregate dimensions.


Abb. 1: Test des linear viskoelastischen Bereichs; Frequenz 1 Hz
Fig. 1: Test of linear region; frequency 1 Hz

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to be built up. Fig. 1 shows the oscillation test as a function of deformation (amplitude) for two different salad dressings. In this diagram, both samples exhibit similar rheological behaviour (pronounced linear-viscoelastic region, indicating a stable system), although the herbs in sample A are in suspension and have settled out in sample B. A further look at the microstructure of the solution, as represented in Fig. 2, shows up the reason for this difference in behaviour. There is no sedimentation if the loss factor $\tan \delta$ is smaller than 1 or if $\delta < 45^\circ$ in the viscoelastic range (see equation 6). The elastic properties dominate and act against the sedimentation force.

2.2 Regeneration

Thixotropic substances in particular show a regenerative phase after their structure has been destroyed, the regeneration behaviour being directly coupled with the flow properties. Regeneration can take seconds or hours, depending on how far the internal structure has been destroyed. To prevent running, application of a coating (e.g. of chocolate) should not take too long. Ketchup and mayonnaise, on the other hand, should quickly revert back to their original firmness soon after application. The regeneration behaviour can be studied using a non-destructive vibration test, as shown in Fig. 3. Here, the structure recovery is detected as a function of time at constant angular frequency ω and deformation γ . Before the actual vibration test, the material being investigated is pre-sheared and the structure destroyed, to simulate actual processing conditions.

Since weak structures are formed in the initial phase, which are immediately destroyed if the deformation is too great, the regeneration behaviour can only be investigated with rheometers capable of producing very small deformations.

2.3 Gelling and curing

Solutions of gelatine and starch in water can occur in states that no longer resemble liquids, but solids. Prepared at elevated temperatures, they become solid on cooling. Pudding retains its shape after cooling and thus has the characteristics of a solid which can be elastically deformed.

During cooling, a coherent network is formed, and the dispersion medium is to be found between the meshes of that network. Fig. 4 shows a typical cooling curve for a gelatine solution, which was cooled from 85 °C down to room temperature. The gelling process starts after about 2.5 minutes and quickly progresses. The

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Fig. 2: Storage and loss moduli of a salad dressing as a function of deformation; frequency 1 Hz

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Fig. 3:
Different regeneration behaviour after preceding shear

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Fig. 4:
Cooling curve for a food gel from $85\text{ }^{\circ}\text{C}$ to room temperature. Frequency 1 Hz, deformation 0.02

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gelling process is complete after 4 minutes and the gel is solid. Margarine exhibits rather different cooling characteristics. Whereas viscous flow is dominant in the case of a gel (see Fig. 4), the elastic properties of this margarine predominate ($G' > G''$), as shown in Fig. 5. Margarine normally consists of 80 % fat which is emulsified in water. This emulsion is stabilised with mono and diacyl glycerides, as well as with crude lethicin. Properties of margarine such as spreadability, a good shelf life and plasticity essentially depend on the type of fat used in its production. Low-fat margarines moreover require a higher emulsifier concentration. The vibration test provides important information as to from which point onward a structure becomes hard, as well as information on processing.

Further information about the internal structure of the substance is provided by dynamic determinations as a function of frequency as shown in Fig. 6. Determinations must be done in the linear-viscoelastic region if meaningful results are to be obtained (see also Fig. 1). The results of such determinations make an important contribution towards clarifying structure [4, 5].

If G' is plotted double-logarithmically against ω , four zones are apparent, which have different characteristics for different substances.

Short-chain linear molecules, which cannot form a network, change direct from the flow zone into the glassy state. The plateau zone indicates a network structure. The stress rate (ω) is faster than the relaxation capacity of the molecules. There is no sliding off of the suspension points and a quasi-solid network is formed. The segments between these suspension points are relatively mobile and when the frequency is increased, the energy is dissipated through the chain segments. The width of this plateau depends upon the molecular weight, and its height on the concentration. The mobility of the chain molecules decreases in the transition zone. The excitation energy is so high that the chain segments can no longer relax completely. The storage

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Fig. 5:
Examination of the crystallisation curve of a commercial grade of margarine, using a vibration test

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
Fig. 6:
Frequency test to characterise the molecular structures

modulus G' becomes bigger and transition into the glassy zone takes place. Here, the chain segments are no longer able to move. The vibration energy is stored through deformation of linkage distances and linkage angles.

3 Conclusions

Foodstuffs show completely different behaviour in the different production stages and applications. Information about thixotropy, rheopexy and viscosity can be obtained from standard measurements.

Dynamic determinations can be used to investigate the behaviour of the various substances in specific applications and to obtain solutions for products and applications. 

- [1] John Dr. Ferry: „Viscoelastics Properties of Polymers“, John Wiley & Sons, Inc., New York, London, Sydney, Toronto (1970)
- [2] W. Flügge: „Viskoelastizität“, Springer Verlag, Berlin, Heidelberg, New York (1975)
- [3] H.A. Barnes: „An Introduction to Rheology“, Elsevier, Amsterdam, Oxford, New York, Tokyo (1989)
- [4] R.B. Bird: „Dynamics of Polymer Liquids, Volume 1: Fluid Mechanics“, 2nd Edition, John Wiley & Sons, New York (1987)
- [5] W. Kulicke: „Fließverhalten von Stoffen und Stoffgemischen“, Hüthig & Wepf Verlag Basel, Heidelberg, New York (1986)